# ILT1000, ILT5000 API Documentation

### **Revision Table**

May 8 2015	v1.5	Initial Release under this form, capturing firmware revisions up to and including 3.0.5.8
May 15 2015	v1.6	Added setwireless command
June 30 2015	v1.7	<ul> <li>Updated getflash to add return value of -513 when asking for a minimum light level that could not be measured in the given range (Rf + voltage sensitivity) due to saturation below the minimum light level</li> <li>added setsampletimetemp functionality to change sample time without saving to flash</li> <li>added usecalfactortemp functionality to change sample time without saving to flash</li> <li>added usefeedbackrestemp functionality to change sample time without saving to flash</li> </ul>
August 11 2015	v1.8	<ul> <li>fixed typo where getambientlevel API was listed as setambientlevel</li> <li>moved setwflisten API to OEM commands</li> </ul>
August 21 2015	v1.9	<ul> <li>Added startpeak, stoppeak, getpeak, startintegrate, stopintegrate, getintegrate</li> <li>Added configbackup, configrestore to OEM commands</li> </ul>

# **ILT1000, ILT5000 API Documentation**

The ILTx000 device uses a simple text-based API for communication with user applications, as well as for use with a basic terminal server.

# **Document Convention**

The API table below contains two columns. In the first column is the API/command name. In the second column is the API/command's definition. The definition contains several elements and convention as follows.

- Actual command and responses will use the Courier New font.
- Command parameters are identified in [brackets], within the "Syntax" description. The brackets themselves are not used in the command syntax, as is illustrated by the "Example command" listed for any command with parameters.
- The "Normal return value(s)" and "Error return value(s)" information is pertinent to the echooff mode of operation (see echoon and echooff commands in the table below). When using echoon mode, return values are self-explanatory.
- An N/A for an "Error return value(s)" indicates that errors are not returned and not expected for this command.
- "Example return" is listed, on some commands, to provide an example of the output. In this field, any text following a dash ("-"), as well as the dash itself, are not part of the actual output and only listed to help describe the output.
- Commands are listed alphabetically.
- The optional "Persist through power-cycle" field will indicate whether or not the configuration is stored in flash memory and "sticks" through power cycles.
- Note that all commands are case sensitive, and always all-lowercase.
- The "Support Starting in Firmware Version" data indicates the firmware version where the command is supported by the manufacturer. The actual command may be present in earlier versions of firmware, but not supported.

# **Product Generations**

The API supports multiple generations of products. At the time of this writing there are three generations: Gen1, Gen2, and Gen3. While there are other differences, the primary changes/additions are as follows:

- Gen2 models add:
  - One additional automatic-gain-control voltage input stage
  - Three programmable feedback resistors (Gen1 was fixed with one value of feedback resistor)
  - A battery back, real-time clock

- o A 4-20mA output
- Gen 3 models add:
  - o An additional, programmable feedback resistor, for a total of four
  - WiFi connectivity (the configuration of which is accomplished with the NetConfig tool)
  - A programmable GPIO pin with current boost, used to control a detector bias voltage (see set0vbias, set5vbias, and getbias)
  - A programmable digital trigger output (see settriggerout)
  - A digital trigger input (see gettriggerin)

# Theory of Operation

Below is the Theory of Operation for interfacing with the device using the API defined further below.

- The device responds to a number of commands (together known as the command-line API) delivered via the USB port.
- The USB port is configured on the device as a USB Serial Port. As a result, interfacing can be done with any standard terminal programs (hyperterminal, putty, MAC terminal window, etc), or direct programming via LabView, .NET, etc.
  - Serial port settings are:
    - Baud Rate: 115200
    - Data Bits: 8
    - Parity bits = None
    - Stop bits = 1
    - Flow Control = None
- The device behaves as follows and in sequence:
  - a. Perform analog/digital conversion and mathematical functions as part of optical level detection
    - i. Buffer up to 4 characters of any incoming commands while performing the processing above
  - b. Check for any characters in the buffer referenced above to indicate an incoming command and process any commands as required. Important programming notes:
    - i. It can take up to 50ms for the system to complete complex analog/digital tasks and process the remainder of the command that was not buffered.
    - ii. The device determines the completion of the command by detecting a "\r" character, i.e. ASCII code 13 decimal.
  - c. Start all over from the top at (a)....
- As a result of the above device behavior, the recommended method to program each device is to:
  - The device operated at 115200 baud, 8 bits, 1 stop bit, no parity bits, no flow control.
  - Always append commands with \r to indicate the command completion.
  - Send the first character of the command, for example the "g" in "getcurrent\r".
  - o Pause 50ms
    - Any setsamplecount settings of more than 50 (2.x firmware) and setsampletime settings of more than 250ms (3.0.5.9 and higher firmware) and include sampling in 50ms or greater chunks of time. During this time, command characters can be dropped. If this is the case, command returns will result in -999, indicating that the command was not understood.
    - This time delay can be lowered when setsampletime and setsamplecount settings are lower.

- Send the remainder of the command, i.e. "etcurrent\r". While there are no pauses required between characters as the device is rapidly draining the incoming serial port at this point, some customers have found more success by inserting a 1ms pause between each character sent within their program.
- Immediately start sensing the response from the device.
  - The device will always send a response to acknowledge the status of the command completion as well as to return values for "get" commands.
  - The response will always be terminated with a "\r\n" (13 decimal, 10 decimal) sequence. This can be used to sense the end of the response and start sending the next command.
    - Regarding time-outs for command responses, "get" commands will typically respond within 100ms. "set" commands that store their configuration in flash memory can take up to 5 seconds to respond. Other special commands like "setfactorydark", "setuserdark", and "captureflash" can take longer to respond.
  - Be careful not to expect an immediate response from the device. For instance, after sending the remainder of the command, i.e. "etcurrent\r", do not immediately check the receive buffer and give-up if data is not sensed. It can take over 1 ms for even the fastest commands to return. It is better (more robust) to (a) send the remainder of the command, (b) continue looking for and processing a response until "\r\n" is sensed, and (c) use a timeout to detect when a response has not returned within some time period.
- Each device is single threaded, meaning that commands and responses need to be processed in sequence. A 2<sup>nd</sup> command cannot, for example, be initiated before the 1st command's response is fully processed. As a result, if a multi-threaded application is accessing the device, a programmatic lock must be placed around device command/response sequences to make sure multiple threads do not attempt to access the device at the same time.
  - If multiple devices are being monitored, a single lock can be used for access to all devices. This is simpler from a coding perspective, but it is not as efficient as it does not allow multiple devices to operate in parallel. For the best performance it is recommended that a per-device lock be established. This allows all devices to be accessed in parallel.
  - As a further performance benefit when monitoring multiple devices, the delay after the first character can be performed in parallel across all devices. For example, if sending "getcurrent\r" to 5 devices, one would:
    - Send "g" to all 5 devices
    - Wait 50ms
    - Send "etcurrent\r" to all devices
    - Process the reply from all devices

Copyright 2015 L&M Instruments LLCPage 4

## API

#### captureflash

Support Starting in Firmware

Version: 3.0.5.8

Support Ending in Firmware

Version: N/A

#### Syntax:

capture flash [trigger type] [trigger level] [minimum light level] [timeout period] [voltage sensitivity] [trigger offset time] [integration time] [save to flash] [save only above minimum light level]

#### What it does:

This command enters a mode where only values are sent back in response to a command. This is the default mode on power up of the device, and the mode used for programming by an application.

#### trigger type:

- 0 = trigger in (start capturing "trigger offset time" milliseconds before trigger in signal)
- 1 = trigger out (start capturing "trigger offset time" milliseconds after trigger out signal)
- 2 = trigger on light level (start capturing after "minimum light level" exceeded)
- 3 = manual trigger (start capturing immediately)

trigger level (applies to trigger type 0 and 1, ignored otherwise):

low = trigger on logic low condition (transition from high to low)

high = trigger on logic high condition (transition from low to high)

#### min light level

Floating point value representing the minimum light level above which capture will start when trigger type is set to 2

#### timeout period

The number of seconds (0 to 300) that the routine will wait for a trigger.

#### voltage sensitivity

- 1 = lowest sensitivity voltage gain (brightest light)
- 2 = medium sensitivity voltage gain
- 3 = highest sensitivity voltage (lowest light level)

#### trigger offset time

The number of milliseconds (0-20) that the trigger will be offset as follows, by trigger type:

- 0 : not applicable
- ${\bf 1}: number\ of\ millise conds\ before\ trigger\ in$
- 2: number of milliseconds after trigger out
- 3: number of milliseconds before minimum light level is sensed

#### integration time

The number of milliseconds (1-40000) that the flash signal will be captured

#### save to flash

0 = do not save to flash (just capture results, see getflash)

1 = save to flash (4096 data points are saved across the integration time, use getlogdata to download)

save only above minimum light level (applicable to trigger type 2)

0 = save all data points

1 = only save data points above the minimum light level

	Ţ
	Normal return value:
	0 on success
captureflash	Error return value(s):
	-500 if missing fields
(Continued)	-501 if bad trigger type. Must be 0-3.
	-502 if bad trigger level. Must be "low" or "high"
	-503 if bad minimum light level. Must be greater than 0.
	-504 if bad timeout period. Must be 0-300.
	-505 if bad voltage sensitivity. Must be 1-3.
	-506 if bad trigger offset. Must be 0-20.
	-507 if bad integration time. Must be 1 – 40000.
	-508 if bad save to flash value. Must be 0 or 1.
	-509 if bad save only above min value. Must be 0 or 1.
	-510 if calibration factor not in use. See usecalfactor.
	-511 if save to flash is set to 1, but data is already in flash. See eraselogdata.
	-512 if timeout period expired before trigger was completed
	-513* if a minimum light level is too high to be measured (device will saturate before
	reaching the level) with given range (Rf + voltage sensitivity)
	*Added in 3.0.6.7
	Persist through power-cycle:
	N/A
	Example command (trigger on light level, "low" trigger level ignored, minimum light level
	= 20e-3 calibrated units, 2 second timeout, Low voltage sensitivity, 5 millisecond trigger
	offset, 40 millisecond integration time after trigger, save data to flash for later retrieval,
	save only the light levels above the minimum trigger level):
	captureflash 2 low 20e-3 2 1 5 40 1 1
clearambientlevel	Syntax:
	clearambientlevel
Support Starting in Firmware	
Version: 3.0.5.8	What it does:
	This command removes any ambient levels configured with setambientlevel.
Support Ending in Firmware	
Version: N/A	Normal return value(s):
	0 on success
	Error return value(s):
	N/A
	Persist through power-cycle:
	No

echooff	Syntax:
	echooff
Support Starting in Firmware	
Version: 2.0.0.3	What it does:
	This command enters a mode where only values are sent back in response to a command.
Support Ending in Firmware Version: N/A	This is the default mode on power up of the device, and the mode used for programming by an application.
	Normal return value:
	0 on success
	Error return value(s):
	N/A
	Persist through power-cycle:
	No
echoon	Syntax:
	echoon
Support Starting in Firmware	
Version: 2.0.0.3	What is does:
	This command enters a verbose mode whereby contextual help is echoed back for each
Support Ending in Firmware Version: N/A	command completion. The mode is useful when interacting with the device from a terminal server or from the CLI program.
	Normal return value:
	0 on success
	Error return value(s):
	N/A
	Persist through power-cycle:
	No

erasecalfactor Syntax: erasecalfactor [calibration number, 1-20] Support Starting in Firmware Version: 2.0.0.3 What it does: This command erases the calibration data associated with the calibration factor. This Support Ending in Firmware includes the calibration factor description, the current-to-irradiance multiplier, and the Version: N/A saturation current. Starting in 3.0.5.3, if the active calibration factor is erased, the cal factor is set to 0 indicating no cal factor in use. Normal return value(s): 0 on success Error return value(s): -500 if missing fields -501 if bad factor number -502 if error erasing flash Persist through power-cycle: Yes Example command: erasecalfactor 5 eraseirrdata Syntax: eraseirrdata Support Starting in Firmware Version: 2.0.0.3 What it does: This command erase any flash memory associated with the tabulated current-to-Support Ending in Firmware irradiance curve established with either setsimpleirrcal or multiple calls to Version: 3.0.5.3 setirrdatapoint. Normal return value(s): 0 on success Error return value(s): N/A Persist through power-cycle: Yes

eraselogdata	Syntax:
	eraselogdata
Support Starting in Firmware	
Version: 2.0.0.3	What it does:
	This command erases any log data that has been stored in the device's flash memory. See
Support Ending in Firmware	<pre>also startlogdata and stoplogdata.</pre>
Version: N/A	
	Normal return value(s):
	0 on success
	Error return value(s):
	-500 if logging is currently active (must use stoplogdata first)
	-501 error erasing the data from flash
	Persist through power-cycle:
	Yes
	163
get100perc	Syntax:
	get100perc
Support Starting in Firmware	goodoopoo
Version: 2.0.0.3	What it does:
	Prior to 3.0.5.3, this command returns the sensor voltage at the time the 100% value was
Support Ending in Firmware	set with set100perc.
Version: N/A	set with section perc.
Version. N//	Starting in 3.0.5.3, this command returns the sensor current at the time the 100% value
	was set with set100perc.
	was set with secrooperc.
	Normal return value(s):
	Prior to 3.0.5.3: Voltage for 100% transmission, volts (firmware >= 2.1.0.0) or microvolts
	(firmware<2.1.0.0)
	Starting in 2.0 E.2. Current for 100% transmission
	Starting in 3.0.5.3: Current for 100% transmission.
	Error return value(s):
	-500 if the 100% value was never set
getambientlevel	Syntax:
	getambientlevel
Support Starting in Firmware	
Version: 3.0.5.8	What it does:
	This command returns the existing ambient light level that is removed from all current
Support Ending in Firmware	and light level readings.
Version: N/A	
,	Normal return value(s):
	Ambient light level, based on detector current
	Error return value(s):
	N/A
	Persist through power-cycle:
	No
	Example return:
	5.981e-6

getambienttemp Syntax: (Gen2, Gen3) getambienttemp Support Starting in Firmware What it does: Version: 2.0.0.3 This command returns the ambient temperature as sensed by the device. This differs from gettemp which returns the internal temperature of the microcontroller. Support Ending in Firmware Version: N/A Normal return value(s): Temperature in degrees F (x100 for firmware < 2.1.0.0). Error return value(s): -500 Command not supported getapiversion Syntax: getapiversion Support Starting in Firmware Version: 2.1.0.0 What it does: This command returns the api version for the firmware. From time-to-time, the CLI API Support Ending in Firmware return values may change the return value or formats. When this happens the API version Version: N/A is increased such that programs can proactively determine the expected return values. API Changes from v1 to v2: get100perc, getvoltage, getvx1, getvx17, getvagc3, getvped, getvref, set100perc all changed from returning microvolts to returning volts in decimal form getcurrent changed from returning picoamps to amps in scientific notation getirradiance changed from returning "pico" values to standard values in scientific notation gettrans changed from returning percent transmission x 10 to percent transmission in decimal form getod changed from returning optical density x 100 to optical density in decimal form getambienttemp changed from returning temp x 100 to tempx1 getlogdata return values changed to reflect unit changes above, with all values returned in scientific notation API Changes from V2 to V3 Added new API's, noted in this document as Supported Starting in Firmware Version 3.0.5.3 or later. set100perc and get100perc changed to return current as opposed to voltage deprecated setclockfreq, setirrdatapoint, storeirrdata, eraseirrdata, setsimpleirrcal, hiddenhelp deprecated setsamplecount (use setsampletime) and getsamplecount (see getsampletime) Overloaded usefeedbackres to return the setting in use if no parameters are provided Normal return value(s): -999 for firmware versions earlier than 2.1.0.0, indicating API version 1. 2 for firmware versions 2.1.0.0 and greater 3 for firmware versions 3.0.5.3 and greater Error return value(s): N/A

getauxserialno	Syntax:
	getauxserialno
Support Starting in Firmware	
Version: 2.0.0.3	What it does:
	This command returns the auxiliary serial number of the device. The serial number is
Support Ending in Firmware	stored in one-time-programmable memory at the time of manufacture.
Version: N/A	
	Normal return value(s):
	Device serial number
	Frank return value(e)
	Error return value(s): -500 device serial number has not been set
	-500 device serial number has not been set
	Example return:
	sn17839-0001
getbias	Syntax:
(Gen3)	getbias
Support Starting in Firmware	What it does:
Version: 3.0.5.4	This command returns the level of bias voltage applied to the detector.
Support Ending in Firmware	Normal return value(s):
Version: N/A	0 indicating no bias is appled
	5 indicating that voltage is applied such that the anode is 5V more negative than the
	cathode
	Error roturn valua(s):
	Error return value(s):
	-501 if not supported on this generation of device

getcalfactor Syntax: getcalfactor [optional, calibration factor, 1-20] Support Starting in Firmware Version: 2.0.0.3 What it does: Without the optional parameter, this command returns the calibration factor currently in Support Ending in Firmware use. With the optional parameter, this command returns the calibration factor details for Version: N/A a particular calibration factor. This includes the following information, separated by space: calibration factor description, current-to-irradiance multiplier(x1000), and saturation current in microamps. A typical use case is to use the command without the optional parameter to determine which calibration factor is in use, followed by issuing the command with the optional calibration factor to determine the details of the calibration factor definition. Normal return value(s): Without the optional command line parameter: 0 if no calibration factor in use 1-20 indicating which calibration factor is in use With the optional command line parameter: [calibration factor description] [current-to-irradiance multiplier(x1000)] [saturation current] Error return value(s): -501 if calibration factor is out of the 1-20 range -502 if the calibration factor is not defined Example command: getcalfactor 1 Example return: calfact1 500 50 getcurrent Syntax: getcurrent Support Starting in Firmware Version: 2.0.0.3 Shortcut: gc (support for this shortcut starting in 3.0.5.4) Support Ending in Firmware Version: N/A What it does: This command returns the sensor current in picoamps. Normal return value(s): Current in amps (firmware>=2.1.0.0) or picoamps (firmware<2.1.0.0) Error return value(s): -500 if there is a voltage saturation, indicating the reading must be discarded

Syntax: getdarkmode getdarkmode Support Starting in Firmware Version: 2.0.0.3 What it does: This command returns the dark mode currently in use by the device. The dark mode can Support Ending in Firmware be NO DARK (there is no consideration for photodiode dark current), FACTORY DARK (the Version: N/A dark current set at the factory), or USER DARK (the dark current set with setuserdark). Return value(s): 0 = NO DARK 1 = FACTORY SET 2 = USER SET Error return value(s): N/A getdatetime Syntax: (Gen2, Gen3) getdatetime Support Starting in Firmware What it does: Version: 2.0.0.3 This command returns the real time clock's date and time in formats as shown below. This format includes a typical date and time format (mm/dd/yyyy hh:mm:ss), followed by Support Ending in Firmware the Epoch time (seconds since 1970). Note that the device is programmed to UTC/GMT Version: N/A time and, as a result, might return a time that does not match the local timezone. Modern computer programming languages will correctly convert the Epoch time data to local time when saving the timestamped log data. This is the case with the Data Logger software. Return value(s): Real time clock time as shown below Error return value(s): -500 Command not supported Example return: 12/05/2013 19:02:05 1386270125 getecal Syntax: (Gen3) getecal Support Starting in Firmware Version: 3.0.5.3 What it does: This command returns the 19 ecal data points. Support Ending in Firmware Version: N/A Normal return value(s): The 19 ecal data points: [1pA val] [5pA val] [10pA val] [50pA [100pA val] [500pA val] [1nA val] [5nA val] [10nA val] [50nA val] [100nA val] [500nA val] [1uA val] [5uA val] [10uA val] [50uA val] [100uA val] [500uA val] [1mA] Error return value(s): N/A

getecaldate (Gen3)	Syntax: getecaldate
Support Starting in Firmware Version: 3.0.5.3  Support Ending in Firmware Version: N/A	What it does:     This command returns the sensor last date electrical cal (ecal) was performed, in epoch time. If an ecal as never performed, this api returns 0.  Normal return value(s):     0 if ecal was never performed     Date in epoch time if ecal was performed  Example return:     1428445793 (indicating April 7 2015 22:29:53 GMT)  Error return value(s):     N/A
getfactorydark  Support Starting in Firmware  Version: 2.0.0.3  Support Ending in Firmware  Version: N/A	Syntax:  getfactorydark  What it does:  This command returns the dark voltage at the various transimpedance amplifier gain stages, as set at the factory.  Normal return value(s):  On Gen1 products: The dark voltage in microvolts, of the two voltage gain stages. On Gen2 products: The dark voltage in microvolts, of the three voltage gain stages for each of three feedback resistor stages On Gen3 products: The dark voltage in microvolts, of the three voltage gain stages for each of four feedback resistor stages  Error return value(s):  -500 if the factory dark voltage has not been set  Example return:  Gen1: 12756 9234 Gen2:
getfeedbackresnumber  Support Starting in Firmware  Version: 3.0.5.8  Support Ending in Firmware  Version: N/A	Syntax:  getfeedbackresnumber  What it does:  This command returns the number, not the value, of the feedback resistor in use.  Normal return value(s):  1 for Gen1 devices (only 1 feedback resistor)  1 – 3 for Gen2 devices  1 – 4 for Gen3 devices  Error return value(s):  N/A

### getflash Syntax: getflash Support Starting in Firmware Version: 3.0.5.8 What it does: This command returns the results of the captureflash command, including the following Support Ending in Firmware Version: N/A Peak: This is the peak light level sensed during the flash capture Average: This is the average light level sensed over the integration period Integration: This is the integration of the light level over the integration period Time Above 10 Percent of Peak: This is the time (in seconds) that the signal was above 10% of the peak value. This can be considered the duration of the pulse. Peak Percent of Range: This is the percentage of the input range (defined by the feedback resistor in use and the voltage sensitivity). A number above 95 would indicate a likely saturation of the signal. The data labels, equal signs, and actual results are all separated by a space, facilitating parsing with a String->Split(' ') command. Normal return value(s): On Gen1 products: The dark voltage in microvolts, of the two voltage gain stages. On Gen2 products: The dark voltage in microvolts, of the three voltage gain stages for each of three feedback resistor stages On Gen3 products: The dark voltage in microvolts, of the three voltage gain stages for each of four feedback resistor stages Error return value(s): The 'Peak' value is returned as -512 if the capture resulted in a timeout before the trigger was sensed. In the response below, the Peak = -5.120e+02 indicates this condition, making the remainder of the values are invalid. Peak = -5.120e+02 Average = 6.168e-06 Integral = 2.527e-07Time-Above-10-Percent-of-Peak = 4.096e-02 Peak-Percent-of-Range = 0Example return: Peak = 1.067e-03 Average = 1.301e-05 Integral = 1.041e-06Time-Above-10-Percent-of-Peak = 7.031e-04 Peak-Percent-of-Range = 97getfeedbackres Syntax: getfeedbackres Support Starting in Firmware Version: 2.0.0.3 What it does: This command returns the value of the feedback resistor, in kilo-Ohms x 10. Support Ending in Firmware Version: N/A Normal return value(s): Transimpedance amplifier feedback resistance in kOhms x 10 Error return value(s): N/A Example return (for a 3K Ohm feedback resistor):

getfriendlyname	I surtou
getifiendfyname	Syntax: getfriendlyname
Support Starting in Firmware	getiliendlyname
Version: 3.0.5.4	What it does:
version. 3.0.3.4	
Cuppert Ending in Firmuras	This command returns the device's "friendly name".
Support Ending in Firmware Version: N/A	
version: N/A	Normal return value(s):
	Device friendly name or "NOT-DEFINED" if the friendly name has not yet been saved
	Error return value(s):
	N/A
getfwversion	Cuntava
getiwversion	Syntax: getfwversion
Support Starting in Firmware	getiwversion
Support Starting in Firmware Version: 2.0.0.3	NA/lash it all a sa
version: 2.0.0.3	What it does:
Constant Fooding in Figure 2	This command returns the firmware version running on the device.
Support Ending in Firmware	
Version: N/A	Normal return value(s):
	Device firmware version
	Error return value(s):
	N/A
	Example return:
	1.3.0.5
getgeneration	Syntax:
geegeneraeron	getgeneration
Support Starting in Firmware	getgeneration
Version: 2.0.0.3	What it does:
Version: 2.0.0.3	
Support Ending in Firmware	This command returns the generation of the device.
	Manusal nations included:
Version: N/A	Normal return value(s):
	-999 for first generation products that did not have this command defined at time of
	release
	1 for first generation products
	2 for second generation products (programmable Rf x 3, realtime clock, ambient temp
	sensor)
	3 for second generation products (programmable Rf x 4, realtime clock, ambient temp
	sensor, wireless)
	Error return value(s):
	N/A
	Funnania matuma
	Example return:
	2
	1

getinfo
Support Starting in Firmware
Version: 3.0.5.3
Support Ending in Firmware
Version: N/A

Syntax:

getinfo

What it does:

This command returns a list, and ever-growing list, of critical device parameters. The command is intended as a quick diagnostic and configuration check, but can be used with parsing software to capture and extract critical configuration information.

Normal return value(s):

See command output for the particular firmware version.

Error return value(s):

N/A

Example return:

```
Base Serial Number = 10002201407300019
Vendor Serial Number = ILT100000002
Model Name = ILT1000-V02
Friendly Name = Right
Generation = 2
FW Version = 3.0.5.8
Dark mode (0=No, 1=Factory, 2=User) = 1
Rf Setting (0=Auto, 1=Rf1, 2=Rf2, 3=Rf3, 4=Rf4) = 0
Rf value R1 (kOhms) = 3
Rf value R2 (kOhms) = 1000
Rf value R3 (kOhms) = 10000
Rf value R4 (kOhms) = 10000000
Rf value in use (KOhms) = 3
eCal: Disabled
eCal values = 1.000e+00 1.000e+00 1.000e+00 1.000e+00
1.000e+00 1.000e+00 1.000e
+00 1.000e+00 1.000e+00 1.000e+00 1.000e+00 1.000e+00
1.000e+00 1.000e+00 1.000e
+00 1.000e+00 1.000e+00 1.000e+00 1.000e+00
ADC Signal Processing: Disabled
Sample time (ms) = 500
Auto Sample Time: Disabled
4-20mA mode (0=Auto/log, 1=min-max/linear, 2=manual) = 0
getvx1 = 0.030337
getvx17 = 0.518344
getvagc3 = 3.300000
Active voltage gain stage = 2
TIA voltage = 0.020656
getcurrent = 6.885e-06
Factory dark = R1 10360 9602 9535 R2 14115 13291 13215 R3
46680 45769 25190
User dark = R1 9735 9607 9564 R2 22885 22746 22670 R3 125018
124804 25190
Logging: Disabled
Wireless Listening: Disabled
```

getintegrate Syntax: getintegrate Support Starting in Firmware Version: 3.0.7.0 What it does: This command returns the integrated light level that started with startintegrate. Support Ending in Firmware The integral is calculated as the average light level since the beginning of Version: N/A startintegrate, multiplied by the time in seconds (resolution of 10ms). Normal return value(s): Integrated light level (units depend on the calibration factor) Error return value(s): -502 if the light level is saturating the detector at the current gain range Example return (for an integrated light level value of 4.712 milli-units): 4.712 - 3getirradiance Syntax: getirradiance Support Starting in Firmware Version: 2.0.0.3 Shortcut: gi (support for this shortcut starting in 3.0.5.4) Support Ending in Firmware Version: N/A What it does: This command returns the irradiance value in user-defined units. The device determines irradiance based on the following priority list: 1. If a calibration factor is in use (see usecalfactor), this current-to-irradiance multiplier value is used. This is the most common method to derive irradiance. 2. If no calibration factor is in use, or the calibration factor has its multiplier value set to 65535 (OxFFFF), the custom calibration table set with either setsimpleirrcal or multiple calls to setirrdatapoint is used. This mode is no longer supported in firmware version 3.0.5.3 and later. Normal return value(s): Irradiance value in user-defined units Error return value(s): -500 if no irradiance calibration data has been set -501 if the detector current is out of range for the current-to-irradiance curve (case #2 above). This value is no longer returned in FW 3.0.5.2 and later. -502 if there is a voltage saturation, indicating the reading must be discarded Example return (for an irradiance value of 7.798 milli-units): 7.798e-3

getirrthresholdlow Syntax: getirrthresholdlow Support Starting in Firmware Version: 2.0.1.0 What it does: This command returns the irradiance value (or light level), below which data will not be Support Ending in Firmware logged. It is set with setirrthresholdlow. Version: N/A Normal return value(s): Minimum irradiance value for data logging, returned in scientific notation. A value of zero (0) can indicate that either the threshold was never set (0 is the default), or it was set to zero. Example return (for an irradiance value of 73.798): 1.25e-4 getlogdata Syntax: getlogdata Support Starting in Firmware Version: 2.0.0.3 What it does: This command returns the log data stored in flash memory. This command can be run Support Ending in Firmware during an active logging session (and datalogging will continue) or after a session is Version: N/A stopped with stoplogdata or, starting in version 3.0.5.8, after captureflash is used with the option to save data to flash. The command will first output three values (total number of values, an integer indicating what information was logged, the logging period, followed by lines of comma-delimited data with date-time-stamp (in seconds since 1970 or "Unix epoch time") followed by all recorded values for that date-time-stamp. Units for Logging Period: Firmware <= 2.0.0.4, with getlogdata used after setlogdata: Logging period in seconds Firmware >= 2.0.0.5, with getlogdata used after setlogdata Logging period in 10ms increments, i.e. 1=10ms, 100=1s Firmware >= 3.0.5.8, with getlogdata used aftercaptureflash Logging period in microseconds, i.e. 10=10us, 1000=1ms Normal return value(s): Total number of Date-Time-Stamp + Value pairs logged Recorded Value Indicator bitmask as follows: 1=Optical Density (x100 for firmware < 2.1.0.0) 2=Percent Transmission (x10 for firmware < 2.1.0.0) 4=Detector Current in amps (picoamps for firmware < 2.1.0.0) 8=Detector Voltage in volts (microvolts for firmware < 2.1.0.0) 16=Device temperature (degrees F) 32=Calibrated Irradiance (see getirradiance) Logging Period (see above regarding units) Seconds Since 1970, value #1, value #2, value #n Seconds Since 1970, value #1, value #2, value #n Seconds Since 1970, value #1, value #2, value #n

Error return value(s):

-500 if no log data present

Notes on returned data values: If the 100% is not set, both Optical Density and Percent Transmission will return Log data does not support negative Optical Density values Calibrated Irradiance will return 0 if there is no calibration data Example return (Notes after the '-' are for documentation purposes and not returned): 5 – total time-stamp + value pairs 4 - detector current - period in seconds (sample every minute) 60 1378738200, 1.595e-9 - 09 Sep 2013 14:50:00 GMT, 1.595 nanoamps 1378738260, 1.346e-9 - 09 Sep 2013 14:51:00 GMT, 1.346 nanoamps 1378738320, 1.456e-9 - 09 Sep 2013 14:52:00 GMT, 1.456 nanoamps 1378738380, 1.748e-9 - 09 Sep 2013 14:53:00 GMT, 1.748 nanoamps 1378738440, 1.637e-9 -09 Sep 2013 14:54:00 GMT, 1.637 nanoamps getmodelname Syntax: getmodelname Support Starting in Firmware Version: 2.0.0.3 What it does: This command returns the model name of the device. Support Ending in Firmware Version: N/A Normal return value(s): Device model name Error return value(s): N/A getod Syntax: getod Support Starting in Firmware Version: 2.0.0.3 What it does: This command returns the optical density, multiplied by 100. Because optical density is a Support Ending in Firmware relative measurement, this command requires that the 100% (0.00 OD) setting is Version: N/A established with set100perc. The maximum optical density returned is based on the device generation as follows: Gen1: 5.000 Gen2: 8.000 Gen3: 11.000 Normal return value(s): Optical density (firmware>=2.1.0.0) or Optical Density x 100 (firmware<2.1.0.0) Error return value(s): -500 if the 100% value has not been previously set with set100perc Example return (for an optical density of 1.070): 1.070

getpeak	Cuntava
getpeak	Syntax: getpeak
Support Starting in Firmware	- Aecheav
Version: 3.0.7.0	What it does:
version. 5.0.7.0	
Cupport Ending in Firmuras	This command returns the peak light level that started when startpeak was executed.
Support Ending in Firmware Version: N/A	
version: N/A	Normal return value(s):
	Peak light level (units depend on the calibration factor)
	Error return value(s):
	-502 if the light level is saturating the detector at the current gain range
	5   1   (5   1   1   1   1   1   1   1   1   1
	Example return (for a peak light level value of 3.465e-4): 3.465e-4
	3.465e-4
getsampletime	Syntax:
accoambic cime	getsampletime
Support Starting in Firmware	geebampieeime
Version: 3.0.5.4	What it does:
VC131011. 3.0.3.4	This command returns the sample time in milliseconds (see setsampletime). Note
Support Ending in Firmware	that if the sample time is set for automatic (set sample time 0) this command will
Version: N/A	return the actual sample in time currently in use.
version. N/A	return the actual sample in time currently in use.
	Normal return value(s):
	Sample time in milliseconds, ranging from 10 to 15000
	Sample time in miniseconds, ranging from 10 to 15000
	Error return value(s):
	N/A
	NA NA
	Example return:
	1000 – indicating 1000ms, or 1 second sample time
	Indicating 1999 in 1999 in the
getserialnumber	Syntax:
	getserialnumber
Support Starting in Firmware	
Version: 2.0.0.3	What it does:
	This command returns the serial number of the device. The serial number is stored in
Support Ending in Firmware	one-time-programmable memory at the time of manufacture.
Version: N/A	, , , , , , , , , , , , , , , , , , ,
-	Normal return value(s):
	Device serial number
	Error return value(s):
	-500 device serial number has not been set
	1 11 11 11 11 11 11 11 11 11 11 11 11 1
	Example return:
	10054201208230245

gettemp	Syntax:
	gettemp
Support Starting in Firmware	
Version: 2.0.0.3	What it does:
Support Ending in Firmware Version: N/A	This command returns the internal temperature of the device microcontroller. Note that, while the detector temperature can be loosely inferred from this, this is not equivalent to the detector temperature. The return value is in degrees F, with any units conversion performed by the application software.
	Normal return value(s):  Device microcontroller temperature, in degrees F
	Error return value(s): N/A
	Example return: 107
gettrans Support Starting in Firmware	Syntax: gettrans
Version: 2.0.0.3	What it does:  This command returns the percent transmission, multiplied by 10. Because percent
Support Ending in Firmware Version: N/A	transmission is a relative measurement, this command requires that the 100% (0.00 OD) setting is established with set100perc.
	Normal return value(s):  Percent Transmission (firmware>=2.1.0.0) or Percent Transmission x 10  (firmware<2.1.0.0)
	Error return value(s): -500 if the 100% value has not been previously set with set100perc
	Example return (for a percent transmission of 67.3): 67.300
gettriggerin (Gen3)	Syntax:  gettriggerin [trigger state] [time out seconds]
Support Starting in Firmware Version: 3.0.5.9	What it does:  This command is for testing of the trigger in line present in some Gen2 and Gen3 devices.
Support Ending in Firmware Version: N/A	The [trigger state] can be "high" (expecting a logic high or "1") or "low" (expecting a logic low or "0"). [time out seconds] is the time to wait for the trigger to appear.
	Normal return value(s):
	0 if trigger not sensed within the time out period 1 if trigger sensed within the time out period
	Error return value(s):
	-500 if missing arguments
	-501 if trigger state is no "low" or "high" -502 if time out period is less than 0 seconds or greater than 300 seconds (5 minutes)

gotugordark	Contract
getuserdark	Syntax:
Company Chambing in Figure 1	getuserdark
Support Starting in Firmware	
Version: 2.0.0.3	What it does:
	This command returns the dark voltage at the various transimpedance amplifier gain
Support Ending in Firmware	stages, as set by the user with setuserdark.
Version: N/A	
	Normal return value(s):
	On Gen1 products: The dark voltage in microvolts, of the two voltage gain stages.
	On Gen2 products: The dark voltage in microvolts, of the three voltage gain stages for
	each of three feedback resistor stages
	_
	On Gen3 products: The dark voltage in microvolts, of the three voltage gain stages for
	each of four feedback resistor stages
	Error return value(s):
	-500 if the user dark voltage has not been set
	-500 II the user dark voltage has not been set
	Example return:
	Gen1:
	13014 9832
	Gen2:
	R1 9735 9607 9564 R2 22885 22746 22670 R3 125018 124804 25190
getvagc3	Syntax:
	getvagc3 (was getvx101 prior to FW 2.0.0.0)
Support Starting in Firmware	
Version: 2.0.0.3	What it does:
	This command returns the voltage at the 3 <sup>rd</sup> automatic gain controller (AGC) stage. For
Support Ending in Firmware	"Gen1" products this is a x101 stage. For "Gen2" products, this is a x131 stage.
Version: N/A	denii producta tina is a xioi stage. For denii producta) tina is a xioi stage.
,	Normal return value(s):
	AGC3 voltage in volts (firmware >= 2.1.0.0) or microvolts (firmware<2.1.0.0)
	AGES Voltage III Volts (IIIIII wate >= 2.1.0.0) of filler ovolts (IIIIII wate < 2.1.0.0)
	Error return value(s):
	N/A
	Example return (1.034 volts):
	1.034054
	1

_	
getvoltage	Syntax: getvoltage
Support Starting in Firmware	your out out of
Version: 2.0.0.3	Shortcut
VEI SIUII. 2.U.U.3	Shortcut:
Cupport Ending in Firm	gv (support for this shortcut starting in 3.0.5.4)
Support Ending in Firmware	
Version: N/A	What it does:
	This command returns the voltage output of the transimpedance amplifier, after it is
	passed through the automatic-gain-control circuit.
	Normal return value(s):
	Detector voltage, in volts (firmware >= 2.1.0.0) or microvolts (firmware<2.1.0.0)
	Error return value(s):
	N/A
	Example return (for 2.415896 volts):
	2.415896
getvoltagestage	Syntax:
6	getvoltagestage
Support Starting in Firmware	
Version: 3.0.5.8	What it does:
6	This command returns the voltage sensitivity stage in use by the device.
Support Ending in Firmware	
Version: N/A	Normal return value(s):
	1 - 3
	Error return value(s):
	N/A
	NYA
getvped	Syntax:
	getvped
Support Starting in Firmware	
Version: 2.0.0.3	What it does:
	This command returns the "pedestal" voltage for the transimpedance amplifier at the x1
Support Ending in Firmware	automatic-gain-control stage. This can be used as a diagnostic to verify this important
Version: N/A	voltage, which should be between roughly 10 and 15 millivolts for Gen1 and Gen2
	products.
	Normal return value(s):
	Normal return value(s):  "Pedestal" voltage in volts (firmware >= 2.1.0.0) or microvolts (firmware<2.1.0.0)
	redestal voltage ili volts (ilililwale >- 2.1.0.0) ol liliciovolts (ilililwale<2.1.0.0)
	Error return value(s):
	N/A
	Example return (11.087 mV):
	0.011087

_	Ţ
getvref	Syntax:
Support Starting in Firmware Version: 2.0.0.3  Support Ending in Firmware Version: N/A	What it does:     This command returns the reference voltage for the device's A/D converter. This command can be used as a diagnostic if there is ever a suspicion that the voltage reference is faulty. This is should be close to, but not typically exactly, 3.3V.  Normal return value(s):     Device A/D converter reference voltage, in volts (firmware >= 2.1.0.0) or microvolts (firmware<2.1.0.0)  Error return value(s):     N/A  Example return:     3.291489
getvx1 Support Starting in Firmware Version: 2.0.0.3 Support Ending in Firmware Version: N/A	Syntax:  getvx1  What it does:  This command returns the voltage at the x1 automatic-gain-control stage.  Normal return value(s):  x1 voltage in volts (firmware >= 2.1.0.0) or microvolts (firmware<2.1.0.0)  Error return value(s):  N/A  Example return (1.543 volts):  1.543087
getvx17 (Gen2, Gen3)  Support Starting in Firmware Version: 2.0.0.3  Support Ending in Firmware Version: N/A	Syntax:  getvx17  What it does:  This command returns the voltage at the x17 automatic-gain-control stage. It only applies to Gen2 and later devices that have this gain stage  Normal return value(s):  x17 voltage in volts (firmware >= 2.1.0.0) or microvolts (firmware<2.1.0.0)  Error return value(s):  -500 if command not supported, i.e. on Gen1 devices  Example return (0.782 volts):  0.782512

help	Syntax:
	help
Support Starting in Firmware	
Version: 2.0.0.3	What it does:
	This command list the most common commands, along with their expected return value.
Support Ending in Firmware	
Version: N/A	Normal return value(s):
	List of common commands and return values.
	Error return value(s):
	N/A
set0vbias	Syntax:
(Gen3)	set0vbias
Support Starting in Firmware	What it does:
Version: 3.0.5.3	This command sets the bias, on the photodiode anode, to zero volts (unbiased).
Support Ending in Firmware	Normal return value(s):
Version: N/A	0 on success
	Error return value(s):
	-500 if not supported, i.e. on Gen1 or Gen2 devices
	Persist through power-cycle:
	No

Cuppert Starting in Firmways	set100perc
Support Starting in Firmware	
Version: 2.0.0.3	What it does:  Prior to 3.0.5.3, this command read the detectors voltage level (and by inference its
Support Ending in Firmware Version: N/A	current and irradiance level) and sets this as the 100% value to use in Optical Density and %Transmission calculations.
	Starting in 3.0.5.3, this command reads the detector current level and sets this as the 100% value to use in Optical Density and %Transmission calculations.
	Normal return value(s):  Prior to 3.0.5.3: The 100% voltage value, in volts (firmware >= 2.1.0.0) or microvolts (firmware<2.1.0.0)
	Starting in 3.0.5.3: The 100% current value, in amperes.
	Error return value(s):
	Prior to 3.0.5.3:
	1 if the value is too low to use as a 100% reference voltage, currently set at < 0.020 V 2 if the value is too high to use as a 100% reference voltage, currently set at > 3.200 V
	Starting in 3.0.5.3: N/A
	Persist through power-cycle: No
	Example return (1.421e-5 current 100% or "full scale"): 1.421e-5
set5vbias (Gen3)	Syntax: set5vbias
Support Starting in Firmware Version: 3.0.5.3	What it does:  This command sets the bias, on the photodiode anode, to -5V volts (negatively biased).
Support Ending in Firmware Version: N/A	Normal return value(s): 0 on success
	Error return value(s):
	-500 if not supported, i.e. on Gen1 or Gen2 devices
	Persist through power-cycle: No

setambientlevel Syntax: setambientlevel Support Starting in Firmware Version: 3.0.5.8 What it does: This command sets the existing current and light-level reading as zero, effectively Support Ending in Firmware eliminating ambient light/power levels from subsequent readings. Version: N/A Normal return value(s): 0 on success Error return value(s): N/A Persist through power-cycle: No setautaveraging Syntax: setautaveraging Support Starting in Firmware Version: 2.0.0.3 What it does: This command sets the amount of averaging performed by the device to be set Support Ending in Firmware automatically based on the transimpedance amplifier voltage as follows: Version: N/A For firmware levels < 2.0.0.0 Voltage > 30mv, averaging equivalent of setlowaveraging 3mv < Voltage < 30mv, averaging equivalent to setmedaveraging Voltage < 3mv, averaging equivalent to sethiaveraging</pre> For firmware levels > 2.0.0.0 Current > 10uA, averaging equivalent of setlowaveraging (0.5 sec) 0.32uA < Current < 10uA, averaging equivalent to setmedaveraging (1.0 Current < 0.32uA, averaging equivalent to sethiaveraging (2.0 sec)</pre> This is the default setting for averaging upon power-up. Normal return value(s): 0 on success Error return value(s): N/A Persist through power-cycle: No for firmware < 3.0.5.3 Yes for firmware >= 3.0.5.3

setcalfactor

Support Starting in Firmware

Version: 2.0.0.3

Support Ending in Firmware

Version: N/A

Syntax:

setcalfactor [calibration factor, 1-20] [desc.] [multiplier x 1000] [saturation current uA]

What it does:

This command defines the particular calibration factor details. Details include the calibration factor description (up to 100 characters), the current-to-irradiance multiplier(x1000), and saturation current in microamps.

If the multiplier value is 0xFFFF (65535) the calibration table as established by setirrdatapoint will be used to convert detector current to irradiance.

Normal return value(s):

0 on success

Error return value(s):

- -500 if missing fields
- -501 if calibration factor is out of the 1-20 range
- -502 if multiplier < 0
- -503 error saving to flash

Persist through power-cycle:

Yes

Example command ("calfact1" description, 0.500 current-to-irradiance factor, 50 uA saturation current):

setcalfactor 1 calfact1 500 50

setcurrentloop (Gen2, Gen3)

Support Starting in Firmware Version: 2.0.0.5

Support Ending in Firmware

Version: N/A

Syntax:

```
setcurrentloop log
setcurrentloop midpoint
setcurrentloop [min picoamps] [max picoamps]
setcurrentloop [0-24]
```

#### What it does:

This command controls the 4-20mA current loop output of devices that support such a current loop.

setcurrentloop log sets the device to output a logarithmic scale current that is in relation to the current sensed by the detector. The transfer function is:

#### Gen2:

4-20mA current = (LOG10(detector current)+8)\*3+5 Detector Current =  $10^{((4-20mA Current)-5)/3-8)}$ 

#### Gen3:

4-20mA current = (LOG10(detector current)+11)\*1+5 Detector Current =  $10^{(([4-20mA Current]-5)/1-11)}$ 

If the detector current is below 10 nA (Gen2) or 1pA (Gen3), the 4-20mA current is set to 4mA. As a result, measurements fluctuating around 10nA (Gen2) or 1pA (Gen3) will result in the 4-20mA output bouncing between 4mA and 5mA. If a calibration factor is set, the 4-20mA output will have a similar relationship to irradiance, but multiplied by the calibration factor.

setcurrentloop midpoint sets the device to output a linear scale where:

4mA = 0 detector current

12mA = the detector current when this command was set

20mA = double the detector current when this command was set

setcurrentloop [min picoamps] [max picoamps] sets the device to output a linear scale where:

4mA = min picoamps 20mA = max picoamps

setcurrentloop [0-24] sets the device to manually output the current indicated, in milliamps.

#### Normal return value(s):

0 on success

#### Error return value(s):

- -500 if command not supported, i.e. on Gen1 devices
- -501 if missing fields
- -502 if there is a bad current loop value (applies to setcurrentloop [0-24])

#### Persist through power-cycle:

Yes for setcurrentloop log, setcurrentloop midpoint, and setcurrentloop [min picoamps] [max picoamps].

No for setcurrentloop [0-24], which is intended as more of a test function.

setdatetime Syntax: (Gen2, Gen3) setdatetime 12/05/2013 19:02:05 Support Starting in Firmware What it does: Version: 2.0.0.3 This command sets the real time clock's date and time. It accepts either of the following formats: Support Ending in Firmware mm/dd/yyyy hh:mm:ss Version: N/A mm/dd/yy hh:mm:ss The device is designed such that this date/time setting is UTC/GMT time. The device stores all date/timestamps in Epoch time format, which is later read out and converted to local time by an application. This is the case with the Data Logger software. Normal return value(s): 0 on success Error return value(s): -500 if command not supported -501 if missing fields -502 if parameters are out of range Persist through power-cycle: Yes, assuming coin/cell battery is in place with adequate charge setfriendlyname Syntax: setfriendlyname [friendly name, up to 30 characters] Support Starting in Firmware Version: 3.0.5.4 What it does: This command sets the "friendly name" of the device. The friendly name is used by Support Ending in Firmware applications to help users easily identify which device is in use. Version: N/A Normal return value(s): 0 on success Error return value(s): -500 if missing fields -501 if error storing value Persist through power-cycle: Yes

sethiaveraging Syntax: sethiaveraging Support Starting in Firmware Version: 2.0.0.3 What it does: This command sets the averaging to be performed by the device to high, resulting in Support Ending in Firmware approximately 1 sample every 2 seconds. This mode provides the highest amount of Version: N/A over- sampling, or noise reduction. This mode is typically used for either detecting very low light levels or for monitoring noisy light sources. Normal return value(s): 0 on success Error return value(s): N/A Persist through power-cycle: setirrdatapoint Syntax: setirrdatapoint [detector current nA] [irradiance x1000] Support Starting in Firmware Version: 2.0.0.3 What it does: This command sets a value pair for a custom detector current-to-irradiance function. Support Ending in Firmware Anywhere from 3 to 100 data points can be set with multiple calls to Version: 3.0.5.4 setirrdatapoint. The device will perform a linear interpolation between data points as needed. This current-to-irradiance function can be used on its own or in conjunction with setcalfactor. Normal return value(s): 0 on success Error return value(s): -500 if missing fields -501 if the maximum data points have already been loaded -502 if a duplicate detector current value already exists in the table Persist through power-cycle: No. See storeirrdata to persist between power cycling. Example command (add a data point for 134 uA = 12.345 irradiance units): setirrdatapoint 134000 12345

setirrthresholdlow	Syntax:
	setirrthresholdlow [calibrated reading]
Support Starting in Firmware	
Version: 2.0.1.0	What it does:
Support Ending in Firmware Version: N/A	This command sets the minimum irradiance value (calibrated reading) to use in conjunction with data logging. When set, and irradiance is being monitored with data logging, data will not be recorded unless the irradiance level meets this threshold.
	Normal return value(s): 0 on success
	Error return value(s):
	-500 if missing fields
	-501 if parameters are out of range:
	Less than 0 or greater than 1 for firmware 2.0.1.0 -> 3.0.5.8
	Less than 0 or greater than 1000 for firmware >= 3.0.5.9
	-502 if error saving to flash memory
	Persist through power-cycle:
	Yes
setlowaveraging	Syntax:
Beerowaveraging	setlowaveraging
Support Starting in Firmware	
Version: 2.0.0.3	What it does:
	This command sets the averaging to be performed by the device to low, resulting in
Support Ending in Firmware	approximately 5 samples per second.
Version: N/A	
	Normal return value(s):
	0 on success
	Ferror return valuadely
	Error return value(s):  N/A
	N/A
	Persist through power-cycle:
	No for firmware < 3.0.5.3
	Yes for firmware >= 3.0.5.3
setmedaveraging	Syntax:
	setmedaveraging
Support Starting in Firmware	
Version: 2.0.0.3	What it does:
Support Ending in Firmware	This command sets the averaging to be performed by the device to medium, resulting in
Version: N/A	approximately 2 samples per second.
version. N/A	Normal return value(s):
	0 on success
	5 5 5466665
	Error return value(s):
	N/A
	Persist through power-cycle:
	No for firmware < 3.0.5.3
	Yes for firmware >= 3.0.5.3

setsamplecount

Support Starting in Firmware

Version: 2.0.0.6

Support Ending in Firmware

Version: 3.0.5.3

Syntax:

setsamplecount [1-200]

What it does:

This command sets the number of samples taken by the Analog-to-Digital converter when reading the voltage from the transimpedance amplifier. This is the immiedate averaging done on the input signal, before additional averaging is done by the <code>set\*averaging</code> commands. Typically a value of 200 is used (the default) for standard operation, or a value of 1 is used for high-speed sampling (up to 100 samples/sec). WARNING: modifying this value will impact the communication between the meter and any device talking to it over the USB port. Talk to the manufacturer for details.

Normal return value(s):

0 on success

Error return value(s):

-500 if missing fields

-501 if value is out of range

Persist through power-cycle:

Yes

Example command (for sampling at 100/sec):

setsamplecount 1

Example command (for default sampling rate):

setsamplecount 200

setsampletime	Syntax:
	setsampletime 0
Support Starting in Firmware	
Version: 3.0.5.4	setsampletime [10 - 15000]
Support Ending in Firmware	What it does:
Version: N/A	This command sets the sample time used by the Analog-to-Digital converter when reading the voltage from the transimpedance amplifier in milliseconds. setsampletime 0 sets the sample time to automatic, resulting in 500ms for high level signals and up to 2 seconds for lower signals.
	For typical monitoring applications, sample time is rarely set below 250ms and typically at 500ms or 1000ms (1s). For high-speed sampling in conjunction with data logging (see startdatalog), sample times as low as 10ms can be used.
	Normal return value(s):
	0 on success
	Error return value(s):
	-500 if missing fields
	-501 if value is out of range
	Persist through power-cycle:
	Yes
	Everyone command (for outs comment times)
	Example command (for auto sample time):  setsampletime 0
	Secsamplecime o
	Example command (for 50ms sample time):
	setsamplecount 50
setsampletimetemp	Syntax:
Cuppert Starting in Firmware	setsampletimetemp 0
Support Starting in Firmware Version: 3.0.6.7	setsampletimetemp[10-15000]
version. 3.0.0.7	2 = c 2 a m b T = c T m e c e m b [ T 0 - T 2000 ]
Support Ending in Firmware	What it does:
Version: N/A	This command performs the same function as setsampletime, but does not persist the setting
-	through power cycles. Because the setting does not need to be saved to flash, along with being
	temporary, it also executes faster.

setsimpleirrcal Syntax: setsimpleirrcal [current-to-irradiance multiplier, x 1,000,000] Support Starting in Firmware Version: 2.0.0.3 What it does: This command sets a simple multiplier to be used to convert detector current to Support Ending in Firmware irradiance. In the absence of factory-defined calibration factors, it is the easiest method Version: 3.0.5.3 to establish a user-defined irradiance value from detector current. The command automatically populates a straight-line (constant) current-to-irradiance relationship using the more advanced setirrdatapoint command. As a result, the storeirrdata command is required to persist this between power cycles. Normal return value(s): 0 on success Error return value(s): -500 if missing fields Persist through power-cycle: No. See storeirrdata to persist between power cycles. Example command (to set a factor of 0.250 to apply to current to derive irradiance): setsimpleirrcal 250000 settriggerout Syntax: (Gen3) settriggerout on settriggerout off Support Starting in Firmware Version: 3.0.5.3 What it does: This command sets the output trigger line to either a logic 1 ("on") or a logic 0 ("off"). This Support Ending in Firmware command is used for testing purposes, with the trigger out circuit being exercised as part Version: N/A of the captureflash command. Normal return value(s): 0 on success Error return value(s): -500 if missing fields -501 if bad arguments Persist through power-cycle: No

setuserdark Syntax: setuserdark Support Starting in Firmware Version: 2.0.0.3 What it does: This command captures the detector voltage signal, at all automatic-gain-control stages, Support Ending in Firmware and stores the values in device flash memory. When used in conjunction with Version: N/A useuserdark, all subsequent readings will have this value removed automatically. Starting in firmware version 2.0.0.5, if the mode is already set to useuserdark then the new values will automatically be applied. Normal return value(s): User dark value, in microvolts, for all gain stages (separated by a space) Error return value(s): -500 if error saving setting to flash memory Persist through power-cycle: Yes Example return: Gen1: 13014 9832 Gen2: R1 9735 9607 9564 R2 22885 22746 22670 R3 125018 124804 25190 setwireless setwireless on (Gen3) setwireless off Support Starting in Firmware Version: 3.0.6.1 What it does: This command enables or disables the wireless functionality. When disabling the Support Ending in Firmware functionality, the system will immediately shut down the wireless part. When enabling Version: N/A the wireless functionality the system requires a power cycle to restore the wireless operation. Normal return value(s): 0 on success Error return value(s): -500 if missing arguments -501 if not supported based on the device generation -502 if bad syntax/arguments Persist through power-cycle: Yes

Support Starting in Firmware

Version: 3.0.7.0

What it does:

This command starts integrating the light level at high-speed to capture fast pulses. The integrated light-level is returned with getintegrate and the integration is halted with stopintegrate. Note that startintegrate always resets the integrated value to 0.

Normal return value(s):
0 on success

Error return value(s):
None

startlogdata

Support Starting in Firmware Version: 2.0.0.3

Support Ending in Firmware

Version: N/A

Syntax:

startlogdata [variable bitmask] [logging period – see below] [seconds since 1970]

#### What it does:

This command defines the logging parameters, and immediately starts logging. This command cannot be run when either a logging session is already active, a session has been stopped (with stoplogdata) but the log data has not yet been erased with eraselogdata, or a captureflash session has completed with saved data that has not yet been erased with eraselogdata. The command takes three parameters, with a space between parameters, as follows:

[variable bit mask]

1=Optical Density (x100)

2=Percent Transmission (x10)

4=Detector Current (picoamps)

8=Detector Voltage (microvolts)

16=Device temperature (degrees F)

32=Calibrated Irradiance (see getirradiance, setirrthresholdlow)

128=Use Real Time Timestamps (as opposed to relative time stamps, Gen2 only)

## For firmware versions 2.0.0.1 and earlier:

[logging period] = desired logging period, in seconds, divided by 10

A 1, for example, would indicate a 10 second delay between logging; a value of 360 would indicate a 3600 second delay, or 1 hour, between log entries. The maximum value is 8640, which is equivalent to 86400 seconds or 1 day. Any values above this will results in a 1 day logging period.

## For firmware versions 2.0.0.2 and later

[logging period] = desired logging period in seconds

A 1, for example, would indicate a 1 second delay between logging; a value of 3600 would indicate a 3600 second delay, or 1 hour, between log entries. The maximum value is 86400, which is equivalent to 86400 seconds or 1 day. Any values above this will results in a 1 day logging period.

# For firmware versions 2.0.1.0 and later

[logging period] = desired logging period, in 10 milliseconds increments.

A 1, for example, would indicate a 10 millisecond delay between logging, or 100Hz logging. A value of 6000 would indicate a 1 minute delay. The maximum value is 8640000, which is equivalent to 86400 seconds or 1 day. Any values above this will results in a 1 day logging period. A value of 0 will default to a 10ms delay. NOTE: there is a known bug that limits the maximum log period to 14.4 minutes. See your product representative for information on an update to resolve this issue if needed.

# [seconds since 1970]

This is "Unix epoch time", with most application coding environments having a mechanism to convert a date-time structure to epoch time. This is ignored, and should be set to 0, when using real-time timestamps.

### Normal return value(s):

0 on success

startlogdata (cont)	Error return value(s):  -500 missing parameters  -501 session already started. Must use stoplogdata and eraselogdata to start new.  -502 errors with the variable bit mask
	Persist through power-cycle: Yes.
	Common usage is to use startlogdata, disconnect the device from the computer, connect the device to a battery or AC-power in the lab or field, and the device will continue logging on power up using time-stamps that are relative to the start time. In this scenario care must be taken to ensure any power down time is negligible to the desired time-stamp accuracy. Gen2 devices have the option of using a battery-back real-time clock (see 128 bitmask value above).
	Example command (to log detector current and device temperature, every 10ms (Firmware 2.0.1.0), starting at 09 Sep 2013 14:50:00 GMT. Note bitmask of 20 is 4/current + 16/temperature):  startlogdata 20 1 1378738200
	Example command. Same as above, but sampling every minute (60 seconds is 6000 10ms increments) and adding 128 to the bit mask for use of a real time clock on a Gen2 device):  startlogdata 148 6000 0
startpeak	Syntax:
Support Starting in Firmware Version: 3.0.7.0 Support Ending in Firmware Version: N/A	what it does:  This command starts peak detection of the light level at high-speed to capture fast pulses of light. The peak light level is returned with getpeak and the peak detection is halted with stoppeak. Note that startpeak always resets the peak value to 0.
	Normal return value(s): 0 on success
	Error return value(s):  None
stopintegrate Support Starting in Firmware	Syntax: stopintegrate
Version: 3.0.7.0	What it does:
Support Ending in Firmware Version: N/A	This command stops the integration started with startintegrate. See also getintegrate.
·	Normal return value(s): 0 on success
	Error return value(s):  -500 if light level integration has not been started with startintegrate

stoplogdata	Syntax:
<u> </u>	stoplogdata
Support Starting in Firmware	
Version: 2.0.0.3	What it does:
	This command stops a log session that was started with startlogdata. Note that,
Support Ending in Firmware	because date-time-stamping is relative to the start time set with startlogdata,
Version: N/A	logging cannot be restarted after a stoplogdata. Instead, log data must first be
	erased using eraselogdata.
	Normal return value(s):
	0 on success
	Error return value(s): -500 if no active logging session to stop
	-300 II IIO active logging session to stop
stoppeak	Syntax:
	stoppeak
Support Starting in Firmware	
Version: 3.0.7.0	What it does:
	This command stops the peak detection started with startpeak. See also getpeak.
Support Ending in Firmware	
Version: N/A	Normal return value(s):
	0 on success
	Error return value(s):
	-500 if peak tracking has not been started with startpeak
	500 ii peak tracking has not been started with 5 car opean
storeirrdata	Syntax:
	storeirrdata
Support Starting in Firmware	
Version: 2.0.0.3	What it does:
	This command stores the custom detector current to irradiance data table, as established
Support Ending in Firmware Version: 3.0.5.3	with multiple calls to setirrdatapoint, in flash memory.
VEISIUII. 3.U.3.3	Normal return value(s):
	0 on success
	5 5.1 5466655
	Error return value(s):
	-500 if no data to store
	Persist through power-cycle:
	Yes

usecalfactor	
usecallactor	Syntax:
Support Starting in Firmware	usecalfactor [calibration number, 0-20]
Version: 2.0.0.3	NATIONAL SE ALONDO
Version, 2.0.0.5	What it does:
Support Ending in Firmware	This command selects a particular calibration factor (see setcalfactor), which in turn will define the detector current-to-irradiance multiplier as well as the detector saturation
Version: N/A	current. A calibration number of 0 is a special case. Using 0 results in no longer using any
	calibration factors. When this is done, getirradiance will only return a value if
	setsimpleirrcal or setirrdatapoint has been used.
	Normal return value(s):
	0 on success
	Error return value(s):
	-500 if missing fields
	-501 if factor number is outside the 0-20 range
	-502 if factor not defined
	-503 if error saving to flash
	Persist through power-cycle:
	Yes (firmware Version 2.0.0.0 and later), No otherwise
	Example command:
	usecalfactor 5
7.5	
usecalfactortemp	Syntax:
Support Starting in Firmware	usecalfactortemp [calibration number, 0-20]
Version: 3.0.6.7	What it does
Version: 3.0.0.7	What it does:
Support Ending in Firmware	This command performs the same function as usecalfactor, but does not persist the setting through power cycles. Because the setting does not need to be saved to flash, along with being
Version: N/A	temporary, it also executes faster.
usefactorydark	Syntax:
	usefactorydark
Support Starting in Firmware	
Version: 2.0.0.3	What it does:
	This command will cause all subsequent readings to remove the factory dark value before
Support Ending in Firmware	presenting any voltage, current, etc. readings. This is the default dark setting upon power
Version: N/A	up.
	Normal return value(s):
	0 on success
	Error return value(s):
	-500 if there has been no factory dark value set
	Persist through power-cycle:
	No

usefeedbackres	Syntax:
(Gen2, Gen3)	usefeedbackres [resistor selection 0-4]
	usefeedbackres
Support Starting in Firmware	
Version: 2.0.0.3	What it does:
	For Gen2 devices, which have 3 available feedback resistors, this command selects the
Support Ending in Firmware	resistor to be used as follows:
Version: N/A	
Tersion 1471	0: have the device automatically select a resistor based on light-level
	1: use feedback resistor #1, usually the lowest value resistor
	2: use feedback resistor #2
	3: use feedback resistor #3
	4: use feedbac resistor #4 (Gen3 only)
	Starting in 3.0.5.4, issuing usefeedbackres without parameters returns the value of the
	feedback resistor in use [0=feedback resistor 1, 1=feedback resistor 2, etc.]
	Starting in 3.0.6.7, issuing usefeedbackres without parameters returns the value of
	usefeedbackres setting [0=auto, 1=fixed to feedback resistor 1, 2=fixed to feedback
	resistor 2, etc.].
	1033001 2, 000.].
	Normal return value(s):
	0 on success (for the first syntax above)
	0-4, depending on the feedback resistor directive (for the second syntax above)
	Error return value(c):
	Error return value(s):
	-500 if missing parameters (pre 3.0.5.4)
	-501 if the device does not support multiple feedback resistors
	-502 it the resistor selection value out of range
	-503 resistor selected, but error saving the change to flash memory
	Persist through power-cycle:
	Yes
usefeedbackrestemp	Syntax:
_	usefeedbackrestemp [resistor selection 0-4]
Support Starting in Firmware	usefeedbackrestemp
Version: 3.0.6.7	
1 0.0.0.0	What it does:
Support Ending in Firmware	This command performs the same function as usefeedbackres, but does not persist the
I I	
Version: N/A	setting through power cycles. Because the setting does not need to be saved to flash, along with
	being temporary, it also executes faster.

usenodark	Syntax:
	usenodark
Support Starting in Firmware	
Version: 2.0.0.3	What it does:
	This command will eliminate any dark current consideration.
Support Ending in Firmware	
Version: N/A	Normal return value(s):
	0 on success
	Error return value(s):
	N/A
	Persist through power-cycle:
	No
useuserdark	Syntax:
	useuserdark
Support Starting in Firmware Version: 2.0.0.3	
version: 2.0.0.3	What it does:
Support Ending in Firmware	This command will cause all subsequent readings to remove the user dark value (see setuserdark) before presenting any voltage, current, etc. readings.
Version: N/A	secuserdark) before presenting any voltage, current, etc. readings.
	Normal return value(s):
	0 on success
	Error return value(s):
	-500 if there has been no user dark value set
	Persist through power-cycle:
	No

# OEM Only

configbackup	Syntax:
	configbackup
Support Starting in Firmware	
Version: 3.0.6.8	What it does:
	This command saves a copy of the device configuration (stored in Sector G) to another
Support Ending in Firmware Version: N/A	area of the device (Sector B). In the even that Sector G becomes corrupted or misconfigured, the factory settings can be restored from Sector B with configrestore.
	Normal return value(s):
	0 on success
	Error return value(s):
	None
	Persist through power-cycle:
	Yes
configrestore	Syntax:
	configrestore
Support Starting in Firmware	
Version: 3.0.6.8	What it does:
Support Ending in Firmware	This command copies the backup configuration (stored in Sector B) to the primary
Support Ending in Firmware Version: N/A	configuration area (Sector G). This is used to restore a device to its factory settings, assuming configbackup was used upon shipment from the factory.
	Normal return value(s):
	0 on success
	Error return value(s):
	None
	Persist through power-cycle:
	Yes

setauxserialno Syntax: setauxserialno [up to 30 non-space characters] Support Starting in Firmware Version: 2.0.0.3 What it does: This command is used by the OEM to set their own serial number. This number is written Support Ending in Firmware into write-once memory on the device, so EXTREME care must be taken to set this Version: N/A properly. Normal return value(s): 0 on success Error return value(s): -500 if a serial number is already programmed -501 if there is an error saving the serial number -502 if the serial number text is missing Persist through power-cycle: Yes Example command (to set the serial number to sn17839-0001): setauxserialno sn17839-0001 setecal Syntax: (Gen3) setecal [1pA val] [5pA val] [10pA val] [50pA val] val] [500pA val] [1nA val] [5nA val] [10nA val] [50nA val] Support Starting in Firmware [100nA val] [500nA val] [1uA val] [5uA val] [10uA val] [50uA [1mA] val] [100uA val] [500uA val] Version: 3.0.5.3 Support Ending in Firmware setecal on Version: N/A setecal off What it does: This first command syntax above inserts calibration points into the device. The procedure involves sourcing the device which current equal to each data point above to determine the value read by the uncalibrated device. The figures are recorded and then entered with the setecal command. The difference between the expected values (above) and the entered values form the basis for the calibration curve. This command syntax also automatically enables ecal. The second command syntax instructs the unit to use the calibrated values. The third command syntax disables the use of the calibrate values. Normal return value(s): 0 on success Error return value(s): -500 if missing fields -501 if the calibration point is too far out of range (+/-4% for expected values >= 50pA, +/-10% on the 5pA and 10pA range). Persist through power-cycle: Yes

setfactorydark

Support Starting in Firmware

Version: 2.0.0.3

Support Ending in Firmware

Version: N/A

Syntax:

setfactorydark

What it does:

This command is used by the manufacturer. It captures the detector signal, at all automatic-gain-control stages, and stores the values in device flash memory. When used in conjunction with usefactorydark, the default dark mode, all subsequent readings will have this value removed automatically. Starting in firmware version 2.0.0.5, if the mode is already set to usefactorydark then the new values will automatically be applied.

Normal return value(s):

Factory dark value, in microvolts, for all gain stages (separated by a space)

Error return value(s):

-500 if error saving setting to flash memory

Persist through power-cycle:

Yes

Example return:

Gen1

11023 10638

Gen2

R1 10360 9602 9535 R2 14115 13291 13215 R3 46680 45769 25190

setfeedbackres Syntax (Gen1 devices): setfeedbackres [feedback resistor value, i.e. Rf, in kOhms x 10] Support Starting in Firmware Syntax (Gen2 devices): Version: 2.0.0.3 setfeedbackres [Resistor 1 kOhms x 10] [Resistor 2 kOhms x 10] [Resistor 3 kOhms x 10] Support Ending in Firmware Version: N/A Syntax (Gen3 devices): setfeedbackres [Resistor 1 kOhms x 10] [Resistor 2 kOhms x 10] [Resistor 3 kOhms x 10] [Resistor 4 kOhms x 10] What it does: This command sets the value of the feedback resistor(s), in kilo-Ohms x 10. On Gen2 devices, which have 3 feedback resistors, this command can take up to three resistor values. On Gen3 devices, which have 4 feedback resistors, this command can take up to three resistor values. Values must be between 10 (1K) and 100000000 (10G). Normal return value(s): 0 on success Error return value(s): -500 if missing fields -501 if error storing value to flash -502 if R1 value is less than 10 Persist through power-cycle: Example Gen1 command (for a 3K Ohm feedback resistor): setfeedbackres 30 Example Gen2 command (for a 3K, 3M, and 10M feedback resistors): setfeedbackres 30 30000 100000 setmodelname Syntax: setmodelname [model name, up to 50 characters] Support Starting in Firmware Version: 2.0.0.3 What it does: This command is used by the manufacturer to set a custom model name. Support Ending in Firmware Version: N/A Normal return value(s): 0 on success Error return value(s): -500 if missing fields -501 if error saving to flash Persist through power-cycle: Yes Example command: setmodelname OpticalMonitor-1200

setwflisten

Support Starting in Firmware

Version: 3.0.5.3

Support Ending in Firmware

Version: N/A

setwflisten on setwflisten off setwflisten tempon setwflisten tempoff

## What it does:

This command disables any listening for commands on the secondary serial port. This is the port that the wireless (WiFly) part uses, but may be used by other devices in the future, including 3<sup>rd</sup>-party connection to the device.

Normal return value(s):

0 on success

Error return value(s):

- -500 if missing arguments
- -502 if bad syntax/arguments

Persist through power-cycle:

Yes for on/off parameters

No for tempon/tempoff parameters